

# INTEX-Ovulationstest<sup>®</sup>

## (OVU-LH)

### Instructions for use

**A rapid, one step test for the qualitative detection of luteinizing hormone (LH) in urine**

#### 1. INTENDED USE

The INTEX-Ovulationstest<sup>®</sup> Test Device (Urine) is a rapid chromatographic immunoassay for the detection of luteinizing hormone (LH) in urine to aid in the early detection of ovulation.

Ovulation is the release of an egg from the ovary. The egg then passes into the fallopian tube where it is ready to be fertilized. In order for pregnancy to occur, the egg must be fertilized by sperm within 24 hours after its release. Immediately prior to ovulation, the body produces a large amount of luteinizing hormone (LH) which triggers the release of a ripened egg from the ovary. This "LH surge" usually takes place in the middle of the menstrual cycle.

The INTEX-Ovulationstest<sup>®</sup> Test Device (Urine) is a complete system to help you predict the time of ovulation, and peak fertility. It is during this fertile time that pregnancy is most likely to occur.

The INTEX-Ovulationstest<sup>®</sup> Test Device (Urine) detects the LH surge in urine, signaling that ovulation is likely to occur in the next 24-48 hours.

Important: The LH surge and ovulation may not occur in all cycles.

#### 2. MATERIALS PROVIDED

- 5 Test devices (each pouch includes 1 Disposable specimen dropper)
- 5 specimen collection containers
- 1 instruction sheet

#### 3. STORAGE AND STABILITY

Store as packaged in the sealed pouch at 4-30°C. The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

#### 4. PRECAUTIONS

- *For in vitro diagnostic use only. Do not use after the expiration date. Do not reuse the test device.*
- *The test device should remain in the sealed pouch until use.*
- *All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.*
- *The test device should be discarded in a proper biohazard container after testing.*
- *The test works only when the test procedures are precisely followed.*

# INTEX<sup>®</sup> Diagnostika

---

- *This test may not be used as a form of birth control.*
- *The test results should not be affected by pain relievers, antibiotics and other common drugs. Medication containing hCG or LH may affect the test and should not be taken while using the INTEX-Ovulationstest<sup>®</sup> Test Device (Urine). In addition, the test will not work properly if you are pregnant, in menopause, or are taking birth control pills.*
- *Keep out of the reach of children.*

## 5. START TESTING

Determine the length of your Menstrual Cycle.

Your Menstrual Cycle length is the number of days from the first day of your period (menstrual bleeding) to the last day before your next period starts. Think back over the last few months to decide what your usual cycle length has been.

Circle your usual cycle length on the WHEN TO START CHART below. Select the number directly underneath. Starting the first day of your last period, count ahead the selected number of days on your calendar. This is the day you should begin testing.

### When to start chart:

Circle your usual cycle length

Cycle length (days)	start with testing (cycle - days)
21	6
22	6
23	7
24	7
25	8
26	9
27	10
28	11
29	12
30	13
31	14
32	15
33	16
34	17
35	18
36	19
37	20
38	21

### Example:

Count ahead this many days STARTING WITH THE FIRST DAY of your LAST period.

Note: If you are unsure about your cycle length, you may want to use your shortest cycle length when reading the chart. If you do this, you may need to test for more than 5 days.

Example: My usual cycle length is 28 days. My last period started on the 6<sup>th</sup> of May. The "When to Start Chart" shows that I will count ahead 11 days beginning with the third. When I count 11 days ahead on the calendar, I find that I will collect and test my urine starting on the 16<sup>th</sup> May". (See "specimen calendar" below).

SU		MO		TU		WE		TH		FR		SA
		1		2		3		4		5		6
7		8		9		10		11		12		13
14		15		16		17		18		19		20
21		22		23		24		25		26		27
28		29		30		31						

Note your personal dates:

Usual cycle length: \_\_\_\_\_

Beginning of the last period at: \_\_\_\_\_

Beginning of the test series (1<sup>st</sup> day) at: \_\_\_\_\_

## 6. SPECIMEN COLLECTION

### Urine assay:

Determine the time you will collect your urine. For best results, collect your urine at about the same time each day, for example, between 10:00AM and 8:00PM. Some women have found that their best specimen is after 12 noon. **Do not collect your first urine after waking up.**

Reduce your liquid intake approximately 2 hours prior to urine collection.

Be sure to write down the date, cycle day and time that you collected your urine in the "Test Results Chart" on the back of this instruction.

You can store your urine for later testing in the day. Your urine can be stored at room temperature for up to 8 hours or in the refrigerator for up to 24 hours. **Do not freeze it.** For best results, test your urine on the same day that it is collected. If you refrigerate your urine, let it reach room temperature before testing (about 30 minutes). Do not shake the container. If a sediment forms at the bottom of the collection container, allow the sediment to settle. Use only urine from the top of the container.

### Specimen storage:

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

### Directions for use:

Allow the test strip, urine specimen and/or controls to equilibrate to room temperature (15-30°C) prior to testing.


Determine the day you will begin testing. (See the above section: "START TESTING").

Bring the pouch to room temperature before opening it. Remove the test device from the sealed pouch and use it as soon as possible.

## 7. TEST PROCEDURE

1

Place the test device on a clean and level surface.




2

Hold the dropper vertically and transfer 5 full drops of urine (approx. 100µl) to the specimen well (S) of the test device, and then start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.

Wait for the red line(s) to appear.

The result should be read at 5 minutes.

It is important that the background is clear before the result is read. Do not interpret the result after 10 minutes.



### Interpretation of results:

**Positive:** Two lines are visible, and the line in test region (T) is darker than the one in the control region (C). It indicates that you will probably ovulate in 24-48 hours.

**Negative:** If two lines are visible, but the line in the test region (T) is lighter than or the same color as the one in the control region (C), or if there is no line in the test region (T). It indicates no LH surge has been detected.

**Invalid:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

### Quality control:

A positive result means your urine contains the level of LH that would indicate you are beginning to ovulate. You can assume that you will ovulate in 24-48 hours.

A negative result means that no LH surge has been detected and you should continue daily testing.

Not every woman ovulates mid-cycle, therefore, you may not see a positive result during the first 5 days of testing. Continue testing with the INTEX-Ovulationstest<sup>®</sup> test Device (Urine).

**RESULTS / QUALITY CONTROL TABLE**  
colour control of test line and control line

DAY	LIGHTER		EQUAL		DUNKLER		DARKER	
	C	T	C	T	C	T	C	T
	Testline T is LIGHTER than the control line C		Testline T is EQUAL (as dark as) the control line C		Testline T is DARKER than the control line C		Testline T is DARKER than the control line C	
	NEGATIVE		POSITIVE		POSITIVE		POSITIVE	
1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
2	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
3	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
4	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
5	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

## 8. QUESTIONS AND ANSWERS

1. *Can I use the INTEX-Ovulationstest Test Device (Urine) to avoid pregnancy?*

No, the test shouldn't be used as a form of birth control.

2. *Do alcohol or common medications affect the test?*

No, but you should consult your physician if you are taking any medication.

3. *How long will the line remain visible?*

The test should be read at 5 minutes for best results. A positive (Surge) result will never disappear. The color line may become darker and a tinted background may appear after several hours. Some negative results may later display a faint second color line because of evaporation from the test region, which prevents complete migration of the test chemicals. Therefore, you should discard the test device once you have read the test.

4. *Once I see a positive result, when is the best time to have intercourse?*

Ovulation is likely to occur within 24-48 hours. This is your most fertile time. Sexual intercourse within this time frame is advised.

5. *I am now using the basal body temperature method (BBT). Does this test replace BBT?*

The shift in basal body temperature primarily indicates that ovulation has already occurred. The INTEX-Ovulationstest Test Device (Urine) indicates that ovulation is about to occur.

6. *I have received a positive result and had intercourse during these fertile days. I have not become pregnant. What shall I do?*

There are many factors that can affect your ability to become pregnant. Often you may need to use the test kit for 3-4 months before achieving pregnancy. You and your partner should consult your physician if pregnancy is not achieved after 3-4 months.

7. *I have had a positive result and had intercourse during these fertile days. I think I may be pregnant. How soon can I find out?*

The INTEX pregnancy products can provide you with results as early as the first day you missed your period.

8. *What time of the day should I perform the test? Do I need to use first morning urine?*

We do not recommend first morning urine, but any other time of day is suitable. For best results, try to collect your urine at approximately the same time each day.






9. *Will the amount of liquid I drink affect the result?*


Heavy intake of fluids prior to testing will dilute the hormone in your urine. We suggest that you limit your fluid intake for about two hours before you collect your urine.

## 9. LITERATURE

1. Bangham, D. R. Acta Endocrinol. 71, 625-637, (1972).
2. Speroff, L., Glass, R. H. Kase N. G. Clinical Gynecologic Endocrinology and infertility, 3<sup>rd</sup> ed., Williams and Wilkins, Baltimore, MD, (1983).
3. France, J. T. In Recent Advances in Obstetrics and Gynaecology Number 14, J. Bonner, ed., Churchill Livingstone, New York, NY, 1982. pp 215-239.
4. Cllins, W. P., Branch, C. M., Collins, P.O., Sallam, H. N., Int J Fert 26, 196-202, (1981).
5. Edwards, R. G. Steptoe, P. C., Purdy, J. M., J Obstet Gynaecol 87,737-756, (1980).
6. Yen, S., Vela, P., and Rankin, J., Journal of Clinical Endocrinology and Metabolism, 30, 435-442, (1970).
7. Engvall, E., Ruoslahti, E. and Engvall, E. J. Immunol. Methods, Vol. 42,11, (1981).

## 10. SYMBOLS

<b>REF</b>	Article number		For single use only
<b>LOT</b>	Lot number		Expiry date
	Storage		Content
<b>IVD</b>	Only for in vitro diagnostics		Instructions for use

<b>INTEX-Ovulationstest<sup>®</sup></b>		<b>REF</b>
	5 test cassettes	OL5505



### INTEX<sup>®</sup> Pharmazeutica AG

Hofackerstrasse 77  
CH-4132 Muttenz  
Tel. (0041) 61 465 90 70  
Fax (0041) 61 465 90 71  
service.ch@intex-diagnostika.com

### INTEX<sup>®</sup> Diagnostika

Hauptstrasse 435  
DE-79576 Weil am Rhein  
Tel. (0049) 7621 940 90 80  
Fax (0049) 7621 940 90 84  
service@intex-diagnostika.com

### INTEX<sup>®</sup> Pharma Hungary Kft.

Árpád u. 3/b. fsz. 1  
HU-1195 Budapest  
Tel. (0036) 1 357 65 98  
Fax (0036) 1 357 65 99  
intex@t-online.hu

[www.intex-diagnostika.com](http://www.intex-diagnostika.com)

Vertrieb Österreich: Mag. Doskar Pharm. Produkte, Schottenring 14, A-1013 Wien, Tel: (0043) 1 535 37 24-0, Fax: (0043) 1 535 37 24-24